

SCITOS – Quick Start Guide

Date: 2013-06-14

Components for your robot

Component	Serial number
SCITOS G5 Base	10388-001
SCITOS G5 HMI	10388-002
Charging Station	10388-003
Laser range finder SICK S300 professional	1320 0191
FLIR Pan-Tilt-Unit PTU D46-17	#007881
Depth camera: Asus Xtion Pro Live	-

EBC Configuration

MainControlUnit EBC:

5V	SCITOS Head
12V	SCITOS Head
24V	SCITOS Head

EBC Card Slot 7:

Port 0	5V	not used
	12V	GigE camera power supply
	24V	not used
Port 1	5V	not used
	12V	FLIR PTU D46-17
	24V	not used

Embedded PC Configuration

CPU:	Intel Core i7-2640M, 2.8GHz
Memory:	4 GB RAM
Harddisc:	sda: 250 GB /dev/sda1 / /dev/sda2 SWAP /dev/sda3 /localhome
Operating system:	Ubuntu 12.04LTS, 32bit, Updates until 2013-05-02
User account:	demo (Password: scitos)
Network configuration:	DHCP on eth0 and eth1
Keyboard settings:	English Keyboard Layout
VNC-Server:	Running at standard port 5900. Password: scitos
CAN adapter:	Integrated MetraLabs-CAN-Adapter (device /dev/ttyUSB2)
CAN bus driver:	<p>If the robot is equipped with a PCAN adapter, the proprietary kernel module <code>pcan</code>. After each kernel update, the module <code>pcan</code> have to be re-compiled:</p> <pre>> sudo -s > cd /usr/local/peak-linux-driver-7.8 > make NET=NO_NETDEV_SUPPORT clean > make NET=NO_NETDEV_SUPPORT > make NET=NO_NETDEV_SUPPORT install</pre> <p>To test the CAN bus interface:</p> <pre>> sudo -s > rmmod pcan > modprobe pcan > cat /proc/pcan → one CAN device (#32) > cat /dev/pcan32 → lots of incoming messages</pre> <p>If your robot has an integrated MetraLabs-CAN-Adapter, you don't have to care about the <code>pcan</code> kernel module.</p>
Touch calibration:	<p>Please use <code>xinput_calibrator</code> to determine the calibration parameters:</p> <pre>> xinput_calibrator</pre> <p>Please copy the new calibration values to the following file:</p> <pre>/usr/share/X11/xorg.conf.d/99-touch</pre>

Software Configuration

Notes:

- The installation of MIRA and MIRA CogniDrive is described in separated document. Please refer to the file `SCITOS-MIRA-InstallGuide-2013-04-30.pdf`.
- At time of the delivery, all software is already installed and ready for use.

MIRA base installation: `/opt/MIRA`

MIRA CogniDrive: `/opt/MIRA-commerical`

MIRA installation repositories:

```
ftp://mira-project.org/repos/  
MIRA-main/ubuntu-1204lts-i686/MIRA-main.repo  
ftp://metralabs-service.com/repos/  
MIRA-commercial/ubuntu-1204lts-i686/MIRA-commercial.repo
```

MIRA license file: `/opt/MIRA-license/*.lic`

Important environment variables in `~demo/.bashrc`:

```
MIRA_PATH=/opt/MIRA:/opt/MIRA-commerical  
LD_LIBRARY_PATH=/opt/MIRA/lib:/opt/MIRA-commerical/lib  
PATH=/opt/MIRA/bin
```

MIRA online documentation:

<http://www.mira-project.org/MIRA-doc/index.html>

Install / update / remove packages using *mirapackage*:

```
> sudo -s  
> source ~demo/.bashrc  
> mirapackage
```

Online documentation of *mirapackage*:

<http://www.mira-project.org/MIRA-doc/MIRAPackagePage.html>

Online documentation of *miracenter*:

<http://www.mira-project.org/MIRA-doc/MIRACenterPage.html>

SCITOS configuration file:

SCITOS robots can be equipped with different kinds of sensors and optional components. For this reason, a global configuration file is used to specify the hardware configuration of the robot.

On your robot, the following configuration file is used:

`/opt/SCITOS/SCITOSRobotAttributes.xml:`

```
<root>
  <var robot="SCITOS-G5" />
  <var frontLaser="SickS300" />
  <var rearLaser="none" />
  <var sonar="false" />
  <var bodyType="normal" />
  <var coverType="2012" />
  <var color="0.10 0.15 0.32" />
  <var hmi="true" />
</root>
```

For more details, please look into the file and read the online documentation:

<http://www.mira-project.org/MIRA-doc/domains/robot/SCITOS/index.html>

Initial SCITOS mapping configuration (please run as user demo):

```
> miracenter -c \
  /opt/MIRA-commerical/domains/robot/SCITOSConfigs/etc/
  SCITOS-mapping.xml
```

- After some initializations, all authorities in the *AuthorityView* should be OK (= green). Otherwise the "Status" should contain a human-readable error message.

Start using CogniDrive:

- Start the SCITOS mapping configuration (see above).
- Record a map of your environment using the *SimpleMapper*.
`http://www.mira-project.org/MIRA-doc/toolboxes/\`
`MapBuilder/index.html`
- Save the new map as "map" in `~demo/mapping/`
- Close the *SimpleMapper* and then close *miracenter*.
- Start the SCITOS navigation configuration:

```
> cd ~demo/mapping
> miracenter -c \
    /opt/MIRA-commerical/domains/robot/SCITOSConfigs/etc/
    SCITOS-Pilot.xml \
    --variables staticMap=/localhome/demo/mapping/map-static.xml
```

Please note, that the path to the map XML files must be an absolute path. Furthermore, this example must be started from the directory, where the maps are located.

- After some initializations, all authorities in the *AuthorityView* should be OK (= green). Otherwise the "Status" should contain a human-readable error message.
- In the view *Channels* you should see a lot of channels, which are updating regularly.
- Open a 3D view: "Window" → "Show view" → "Visualization – 3D"
- Put the following channels into the 3D view (using drag & drop from the *Channels* view):
`/maps/static/Map`
`/robot/Laser`
`/robot/navigation/Path`
- Go to the "Visualization Control" window a using the "+" Button to add a "Model – Rigid Model" visualization.
- Now you should see the robot, the map and the laser range finder data in the 3D view. Use your mouse to move around in the 3D view ("Orbit mode"):
Left Button → Rotation
Right Button → Moving
Scroll Wheel → Zooming
- The 3D view contains two buttons "Localize" and "Goal" to localize the robot and to set goal position. See the tool tips of this tools for some instructions. More details can be found in the online documentation of *miracenter*. Before using this tools, please record a map of your environment!
- Open the *DriveView*: "Window" → "Show view" → "Robot – DriveView". With this view, you should be able to move the robot round using the keyboard. The *DriveView* also contains button for reset the bumper ("Reset Motorstop") and to enable/disable the *FreeRunMode* of the robot.